



Vcheck
cPL



A Single-Step Testing Procedure for Pancreatitis in Dogs

Acute pancreatitis in dogs can be extremely painful. It can become life-threatening by affecting multiple organs, if the inflammation spreads. Though common in dogs, it is often difficult to diagnose. Vcheck cPL has the most sensitive and specific test for canine pancreatitis on the market. It is accurate and easy to follow, minimizing the chance of user error and improving the reproducibility of results.

Clinical Applications

- To diagnose acute pancreatitis when nonspecific symptoms occur
- Monitoring response to therapy by serial checking for, and evaluating, treatment efficacy
- To assess secondary damage to the pancreas

Specifications

Species	Canine
Sample Type	Serum 25 μ l
Measurement	Quantitative
Range	50 - 2,000 ng/ml
Testing Time	5 minutes
Storage Condition	1 - 30° C

Simple Testing Procedure



Dilute Sample

Use a 25 μ l pipette to draw 25 μ l of serum and add to the assay diluent tube.



Mix

Use a 100 μ l pipette to mix the sample with diluent by pipetting 5 - 6 times.




Measure

Add 100 μ l of the mixed sample to the sample well of the test device and press [START].

Product Name	Product Number	Product Type	Packing Unit
Vcheck cPL 2.0	VCF129DD	Device	10 Tests/Kit

A Closer Look: cPL

Vcheck cPL has a high correlation with company “I” laboratories. This analyzer allows in-clinic, quantitative performance measurements for an accurate diagnosis of pancreatitis.



Specific Clinical Application

With the onset of acute pancreatitis, timely and accurate testing drastically improves the likelihood of proper treatment. Time is critical when analyzing and treating a canine in this situation. The Vcheck cPL analyzer provides timely analysis by providing rapid, in-clinic testing, with reproducible and accurate results.

PLI Testing

< 200 ng/ml

Pancreatitis very unlikely

200-400 ng/ml

Equivocal

> 400 ng/ml

Pancreatitis

*Internal Evaluation Data



For More Information on
Vcheck V200 or V2400
analyzers visit:
bionote.com
customerservice@bionote.com
800-727-5169

LEARN MORE
about the cPL test

