

Summary of Bionote Study: Assessment of Vcheck[®] analyzer for rapid progesterone concentration measurement including recommendations for achieving the optimal breeding time

Supphathat Wutthiwitthayaphon,¹ Thanikran Suwannachote,² Saengtawan Arayatham,^{3,4} Wisut Prasitsuwan,³ and Sakchai Ruenphet^{1,3}

1. Department Immunology and Virology, Veterinary Medicine Faculty, Mahanakorn University of Technology, 140 Cheum-Sampan Rd. Nong Chock, Bangkok, Thailand
2. Clinic for Small Domestic Animals and Radiology, Veterinary Medicine Faculty, Mahanakorn University of Technology, 140 Cheum-Sampan Rd. Nong Chock, Bangkok, Thailand
3. Master of Science Program in Animal Biotechnology, Veterinary Medicine Faculty, Mahanakorn University of Technology, 140 Cheum-Sampan Rd. Nong Chock, Bangkok, Thailand
4. Clinic of Obstetrics, Gynecology and Animal Reproduction, Veterinary Medicine Faculty, Mahanakorn University of Technology, 140 Cheum-Sampan Rd. Nong Chock, Bangkok, Thailand

Background

Serum progesterone concentration plays critical role in determining the optimal breeding time in bitches and diagnosing reproductive-related issues. This study aimed to conduct a comparative analysis of serum progesterone results obtained from commercial point-of-care immunological analyzers, namely, Vcheck[®], with those obtained using chemiluminescent microparticle immunoassay (CMIA). Our overarching goal was to evaluate these analyzers' accuracy and establish standardized guidelines for optimal breeding timing.

Introduction

The optimal breeding time in bitches requires an assessment of serum progesterone concentrations^[1]. However, this evaluation serves multiple essential purposes, including identifying reproductive irregularities, such as hypoluteoidism^[2], and confirming luteolysis before parturition^[3-5]. Notably, a characteristic increase in serum progesterone concentrations during the estrus period often exceeding 1 ng/mL, indicates significant hormonal changes. In bitches, ovulation typically occurs 36–50 h after the luteinizing hormone (LH) peak^[6], which correlates with serum progesterone concentrations of approximately 2.02–0.18 ng/mL^[7]. These concentrations then escalate to a range of 4.00–10.00 ng/mL on the day of ovulation^[8], indicating a significant hormonal shift and the onset of ovulation. Intriguingly, despite this range, Seefeldt *et al.*^[9], Marseloo *et al.*^[10], and Mir *et al.*^[11] suggested a serum progesterone concentration of 5.00–8.00 ng/mL, introducing contrasting perspectives on determining the ovulation.

Various techniques, such as radioimmunoassay (RIA)^[12,13], liquid chromatography-tandem mass spectrometry (LC-MS)^[14,15], and chemiluminescence immunoassay (CLIA)^[14,16], are used to measure serum progesterone concentration in veterinary medicine. An accurate and reliable alternative is the CLIA method, adept at assaying serial blood samples with guaranteed safety, speed, accuracy and repeatability^[17]. CLIA addresses the drawbacks associated with RIA and enzyme immunosorbent assay and serves as a robust solution for precise serum progesterone monitoring, aiding in accurate ovulation prediction and confirmation. Moreover, recent advances, such as point-of-care analyzers, have transformed the measurement of progesterone. These advancements, including rapid fluorescence immunochromatography assay, surface plasmon field-

enhanced fluorescence spectroscopy, lateral flow immunochromatography and competitive enzyme-linked fluorescence assay^[18-21], enhance serum progesterone monitoring, ultimately improving veterinary practice in managing reproductive processes.

However, differences in serum progesterone levels are due to the use of different laboratory techniques and variations between bitches. As a result, the precise identification of the ideal mating period necessitates the collection of multiple consecutive blood samples during both the proestrus and estrus phases, which can then be compared with established gold standards or reference laboratory procedures. This investigation involves a comparative analysis of serum progesterone findings derived from a commercial point-of-care analyzer, namely Vcheck[®], in contrast to those acquired through chemiluminescent microparticle immunoassay (CMIA) using the same serum samples.

Materials and Methods

Ninety-four serum samples from bitches were analyzed using the Vcheck[®] analyzer and compared with CMIA. Thorough documentation included the mean, standard deviation, 95% confidence interval (CI), and minimum and maximum values of serum progesterone concentrations. Furthermore, Pearson's correlation coefficient, Lin's concordance correlation coefficient, and the bias correction factor were meticulously recorded.

Results

The mean progesterone concentration measured using the Vcheck[®] analyzer was significantly lower than that measured using CMIA, with a mean difference of 1.26 ng/mL of serum. The Bias correction factor was 0.935, which was nearly 1.00, indicating that the line of best-fit was on the perfect line of agreement, providing insight into the measurement accuracy. Pearson's correlation coefficient, a measure of precision, was also close to 1 (0.939), confirming the reliability of the data. Furthermore, Lin's concordance correlation coefficient was 0.877, indicating a fair overall agreement between the Vcheck[®] and CMIA methods. These results support the validity of the Vcheck[®] analyzer's results. The present study was developed by aligning with established CMIA guidelines and adapting them using the range and 95% CI derived from each set of results, ensuring a standardized and rigorous approach.

Table 1: Provides a comprehensive overview of the means, SD, 95% CI, and range of serum progesterone concentrations. These values were determined through measurements taken during various phases, including early proestrus, LH peak, pre-ovulation, ovulation, post-ovulation and all bitch periods. This detailed presentation enables us to understand the distribution of progesterone levels in different phases.

Mean, SD, 95% CI, Min, and Max value for serum progesterone concentration with quantification using the CMIA and Vcheck® for estimates during the early proestrus, LH peak, pre-ovulation, ovulation, post-ovulation and all periods of the bitch.

Period	CMIA			Vcheck®		
	Mean ± SD	95% CI	Min-Max	Mean ± SD	95% CI	Min-Max
Proestrus	1.24 ± 0.37	1.11 - 1.38	1.00 - 1.98	1.36 ± 0.54	1.16 - 1.56	1.00 - 2.81
LH peak	2.65 ± 0.36 ^a	2.35 - 2.95	2.17 - 2.98	1.88 ± 0.55 ^b	1.42 - 2.33	1.00 - 2.73
Pre-ovulation	3.71 ± 0.52	3.31 - 4.10	3.00 - 4.44	3.35 ± 1.53	2.17 - 4.53	1.77 - 7.01
Ovulation	6.79 ± 1.32 ^a	6.31 - 7.27	5.09 - 9.54	5.36 ± 2.48 ^b	4.47 - 6.25	1.00 - 10.60
Post-ovulation	16.47 ± 7.18 ^a	12.32 - 20.61	10.11 - 30.00	11.7 ± 5.90 ^b	10.26 - 14.85	5.77 - 21.27
All period	5.75 ± 5.79 ^a	4.57 - 6.94	1.00 - 30.00	4.50 ± 4.39 ^b	3.60 - 5.40	1.00 - 21.27


^{a,b}values within the same row with different superscripts mean statistically significant difference (p < 0.05). SD=Standard deviation, CI=Confidence interval, Min=minimum, Max=maximum, CMIA=Chemiluminescence microparticle immunoassay

We found an intriguing finding when we compared these values: The average values of all samples obtained using Vcheck® were significantly lower than those obtained using CMIA, with an average difference of 1.26 ng/mL (Table-2 and Figure-1). This difference is not only statistically significant but also significantly impacts the accuracy and precision of the measurements.

Table 2: Measure of agreement: Concordance correlation coefficient, Pearson's' correlation coefficient, and bias correction factor.

95% limits of agreement (Bland and Altman)			Lin's concordance correlation coefficient	95% confidence interval		Pearsons' correlation coefficient	Bias correction factor
Average difference	Lower	Upper		Lower	Upper		
1.26	-3.16	5.68	0.877	0.834	0.910	0.939	0.935

Figure 1: Passing-Bablok regression plot depicting serum progesterone measurements by Vcheck® against the chemiluminescence microparticle immunoassay serum progesterone measurements. The thin line is the identity line. The thick line represents the regression line and the dotted lines represent its 95% confidence interval.



Conclusion

The Vcheck® analyzer provides a rapid assessment of serum progesterone concentration in bitches, with results comparable to those measured using the CMIA technique. However, when considering the use of the Vcheck® analyzer, it is recommended that the results should be interpreted carefully and the interpretation guidelines should be followed. In conclusion, Vcheck® provides a reliable and convenient method for veterinarian practitioners to measure canine progesterone levels in a clinical/hospital setting.

Sources

1. Conley A.J, Gonzales K, Erb H.N, Christensen B.W. Progesterone analysis in canine breeding management. *Vet. Clin. North Am. Small Anim. Pract.* 2023;53(5):931–949. [PubMed] [Google Scholar]
2. Tibold A, Thuroczy J. Progesterone, oestradiol, FSH and LH concentrations in serum of progesterone-treated pregnant bitches with suspected luteal insufficiency. *Reprod. Domest. Anim.* 2009;44(2):129–132. [PubMed] [Google Scholar]
3. De Cramer K.G.M, Nöthling J.O. The precision of predicting the time of onset of parturition in the bitch using the level of progesterone in plasma during the preparturient period. *Theriogenology.* 2018;107(11):21–218. [PubMed] [Google Scholar]
4. Roos J, Maenhoudt C, Zilberstein L, Mir F, Borges P, Further E, Niewiadomska Z, Nudelman N, Fontbonne A. Neonatal puppy survival after planned caesarean section in the bitch using aglepristone as a primer: A retrospective study on 74 cases. *Reprod. Domest. Anim.* 2018;53(suppl 3):85–95. [PubMed] [Google Scholar]
5. Nöthling J.O, Jooné C.J, Hegarty E, Schooley E.K, De Cramer K.G.M. Use of a point-of-care progesterone assay to predict onset of parturition in the bitch. *Front. Vet. Sci.* 2022;9:914659. [PMC free article] [PubMed] [Google Scholar]
6. Concannon P.W. Canine physiology of reproduction. In: Burk J, editor. *Small Animal Reproduction and Infertility.* Philadelphia USA: Lea and Febiger; 1986. pp. 23–77. [Google Scholar]
7. Kutzler M.A, Mohammed H.O, Lamb S.V, Meyers-Wallen V.N. Accuracy of canine parturition date prediction from the initial rise in preovulatory progesterone concentration. *Theriogenology.* 2003;60(6):1187–1196. [PubMed] [Google Scholar]
8. Edens M.S.D, Heath A.M. Breeding management in the bitch and queen. In: Root Kustritz M.V, Messonnier S.P, editors. *Small Animal Theriogenology.* Vol. 2. Saint Louis, USA: Butterworth-Heinemann; 2003. p. 33. [Google Scholar]
9. Seefeldt A, Schone J, Brussow N, Bunck C, Hoppen H.O, Beyerbach M, Günzel-Apel A.R. Relevance and accuracy of ovulation timing with regard to prediction of parturition in the dog. *Tierarztl. Prax. K. H.* 2007;35(3):188–192. [Google Scholar]
10. Marseloo N, Fontbonne A, Bassu G, Rivire S, Leblanc B, Rault D, Biourge V, Chastant-Maillard S. Proceedings of the 5th International Symposium on Canine and Feline Reproduction, Embu´ das Artes, SP, Fontbonne, Brazil. 2004. Comparison of Ovarian Ultrasonography with Hormonal Parameters for Determination of the Time of Ovulation in Bitches; pp. 75–77. [Google Scholar]
11. Mir F, Billault C, Fontaine E, Sendra J, Fontbonne A. Estimated pregnancy length from ovulation to parturition in the bitch and its influencing factors: A retrospective study in 162 pregnancies. *Reprod. Domest. Anim.* 2011;46(6):994–998. [PubMed] [Google Scholar]
12. Skenandore C.S, Pineda A, Bahr J.M, Newell-Fugate A.E, Cardoso F.C. Evaluation of a commercially available radioimmunoassay and enzyme immunoassay for the analysis of progesterone and estradiol and the comparison of two extraction efficiency methods. *Domest. Anim. Endocrinol.* 2017;60:61–66. [PubMed] [Google Scholar]
13. Tal S, Mazaki-Tovi M, Druker S, Novak S, Raz T, Aroch I. Evaluation of two chemiluminescent assays compared with radioimmunoassay for serum progesterone measurement in bitches. *Theriogenology.* 2020;147:116–123. [PubMed] [Google Scholar]
14. Trabert B, Falk R.T, Stanczyk F.Z, McGlynn K.A, Brinton L.A, Xu X. Reproducibility of an assay to measure serum progesterone metabolites that may be related to breast cancer risk using liquid chromatography-tandem mass spectrometry. *Horm. Mol. Biol. Clin. Investig.* 2015;23(3):79–84. [PMC free article] [PubMed] [Google Scholar]
15. Sasaki M, Ochiai H, Takahashi K, Suzuki R, Minato K, Fujikata A. Development and validation of LC-MS/MS assay for the quantification of progesterone in rat plasma and its application to pharmacokinetic studies. *Drug Res. (Stuttg)* 2015;65(9):484–489. [PubMed] [Google Scholar]
16. Nöthling J.O, De Cramer K.G.M. Comparing the values of progesterone in the blood of bitches as measured with a chemiluminescence immunoassay and a radioimmunoassay. *Reprod. Domest. Anim.* 2018;53(5):1136–1141. [PubMed] [Google Scholar]
17. Boland P, Potts G, Chapwanya A, Vangen M. The Whelping Rate and Litter Size of 157 Surgical Artificial Inseminations with Frozen Semen in Racing and Coursing Greyhounds. *ATTI IV EVSSAR Congress.* 2004:37–38. [Google Scholar]
18. Kunanusont N, Punyadarsaniya D, Ruenphet S. Accuracy and precision guidelines for optimal breeding time in bitches using in-house progesterone measurement compared with chemiluminescent microparticle immunoassay. *Vet. World.* 2021;14(3):585–588. [PMC free article] [PubMed] [Google Scholar]
19. Fontbonne A, Maenhoudt C, Thoumire S, Roos J, Niewiadomska Z, Robiteau G, Rousseliere E, Buronfosse T. Evaluation of surface plasmon field-enhanced fluorescence spectroscopy for rapid measurement of progesterone concentration in bitches. *Am. J. Vet. Res.* 2021;82(5):417–424. [PubMed] [Google Scholar]
20. Milani C, Boscato E, Gabai G, Badon T, Schrank M, Sontas H, Romagnoli S, Mollo A. Analytical and clinical performance of a fluorescence enzyme immunoassay for progesterone and determination of ovulation day in bitches. *J. Vet. Diag. Invest.* 2022;34(6):977–982. [PMC free article] [PubMed] [Google Scholar]
21. Hussein H.A, Sculer G, Wehrend A. Comparison of three progesterone quantification methods using blood samples drawn from bitches during the periovulatory phase. *Vet. World.* 2022;15(1):119–123. [PMC free article] [PubMed] [Google Scholar]
22. Luz M.R, Bertan C.M, Binelli M, Lopes M.D. Plasma concentrations of 13,14-dihydro-15-keto prostaglandin F₂-alpha (PGFM), progesterone and estradiol in pregnant and nonpregnant diestrus cross-bred bitches. *Theriogenology.* 2006;66(6–7):1436–1441. [PubMed] [Google Scholar]
23. Okkens A.C, Teunissen J.M, Van Osch W, van den Brom W.E, Dieleman S.J, Kooistra H.S. Influence of litter size and breed on the duration of gestation in dogs. *J. Reprod. Fertil. Suppl.* 2001;57:193–197. [PubMed] [Google Scholar]
24. Steckler D, Nöthling J.O, Harper C. Prediction of the optimal time for insemination using frozen-thawed semen in a multi-sire insemination trial in bitches. *Anim. Reprod. Sci.* 2013;142(3–4):191–197. [PubMed] [Google Scholar]
25. Rota A, Charles C, Cucuzza A.S, Pregel P. Diagnostic efficacy of a single progesterone determination to assess full-term pregnancy in the bitch. *Reprod. Domest. Anim.* 2015;50(6):1028–1031. [PubMed] [Google Scholar]
26. Schmicke M, Urhausen C, Wolf K, Schmidt S, Günzel-Apel A.R. Evaluation of the blood progesterone concentration in the bitch measured by chemiluminescence immunoassay at the day of ovulation. *Tierarztl. Prax. Ausg. K Kleintiere Heimtiere.* 2016;44(5):317–322. [PubMed] [Google Scholar]
27. Volkmann D.H. The effects of storage time and temperature and anticoagulant on laboratory measurements of canine blood progesterone concentrations. *Theriogenology.* 2006;66(6–7):1583–1586. [PubMed] [Google Scholar]
28. Bolelli G, Muti P, Micheli A, Sciajno R, Franceschetti F, Krogh V, Pisani P, Berrino F. Validity for epidemiological studies of long-term cryoconservation of steroid and protein hormones in serum and plasma. *Cancer Epidemiol. Biomarkers Prev.* 1995;4(5):509–513. [PubMed] [Google Scholar]
29. Key T.J, Pike M.C, Baron J.A, Moore J.W, Wang D.Y, Thomas B.S, Bulbrook R.D. Cigarette smoking and steroid hormones in women. *J. Steroid Biochem. Mol. Biol.* 1991;39(4A):529–534. [PubMed] [Google Scholar]
30. Zuercher J, Boes K.M, Balogh O, Helms A.B, Cecere J.T. Comparison of a point-of-care analyzer with a chemiluminescent immunoassay for serum progesterone measurement in breeding management of the bitch. *Front. Vet. Sci.* 2021;8:660923. [PMC free article] [PubMed] [Google Scholar]
31. Taylor R. Interpretation of the correlation coefficient: A basic review. *J. Diagn. Med. Sonogr.* 1990;6(1):35–39. [Google Scholar]

VIEW THE FULL
STUDY HERE

