



Vcheck V200



Your Best Treatment Plan is in the Numbers

Affordable, rapid and proven, Vcheck is the only in-clinic immunofluorescence assay analyzer that captures 16 biomarkers, including Canine NT-proBNP. Vcheck improves patient care by delivering actual quantitative results in minutes, eliminating estimates, ambiguity and the need for reference lab confirmation. As a multiple-test analyzer, Vcheck can be used in-clinic as a complement to your point-of-care lab for time-saving diagnosis and faster initiation of treatment plans.

Precise quantitative result

No positives, negatives or estimates

Significant cost savings per test

The most affordable analyzer available

Easy user interface

Optical Character Recognition (OCR) reads handwriting directly from test

Auto-coding system with 2D barcode technology

Randomly accessible without any pre-procedure

Clinic-friendly size

Specifications

Model	Vcheck V200
Test capacity	1 test at a time
Power	AC/DC adaptor 120V
Display	7" Color Touch Screen
Printer	Built-in
Connectivity	HL7 v2.6(PCD-01) / POCT1-A
Dimension	7.8" Wide x 9.4" Deep x 8" High
Weight	5.5 lbs

Simple Testing Procedure



1. Dilute the sample



2. Add diluted sample into a test device



3. Check the test result

Available Test Reagents

Cardiac

Detection and monitoring of heart disease in dogs and cats, including MMVD and DCM

- Canine NT-proBNP
- Feline NT-proBNP
- Canine Tnl
- Feline Tnl

Acute Inflammation

Earlier detection of acute inflammation in dogs, cats and equine

- Canine CRP
- Feline SAA
- Equine SAA

Pancreatitis

Definite diagnosis of Pancreatitis in dogs and cats

- cPL (canine Pancreas-specific Lipase)
- fPL (feline Pancreas-specific Lipase)

Coagulation

Highly sensitive marker for thromboembolism

- Canine D-dimer

Hormones

A hormone marker for Cushing's and Addison's diseases

- Canine Cortisol

A hormone marker for diagnosis of thyroid gland diseases

- T4 (Thyroxine)
- Canine TSH

For determining optimal breeding dates

- Canine Progesterone
- Equine Progesterone

Antibody

Assessment and monitoring of immune levels of neonatal foals.

- Foal IgG

Product Name	Product Number	Product Type	Packing Unit
Vcheck Canine NT-pro BNP	VCF132DC	Device	5 Tests/Kit
Vcheck Feline NT-pro BNP	VCF130DC	Device	5 Tests/Kit
Vcheck Canine CRP 2.0	VCF109DD	Device	10 Tests/Kit
Vcheck Feline SAA 3.0	VCF138DD	Device	10 Tests/Kit
Vcheck cPL 2.0	VCF129DD	Device	10 Tests/Kit
Vcheck fPL 2.0	VCF127DD	Device	10 Tests/Kit
Vcheck D-dimer	VCF107DD	Device	10 Tests/Kit
Vcheck cCortisol	VCF105DD	Device	10 Tests/Kit
Vcheck T4	VCF106DD	Device	10 Tests/Kit
Vcheck Canine TSH	VCF118DC	Device	5 Tests/Kit
Vcheck cProgesterone	VCF122DD	Device	10 Tests/Kit
Vcheck Canine Tnl	VCF137DC	Device	5 Tests/Kit
Vcheck Feline Tnl	VCF139DC	Device	5 Tests/Kit
Vcheck Equine SAA	VCF141DC	Device	10 Tests/Kit
Vcheck eProgesterone	VCF142DC	Device	5 Tests/Kit
Vcheck Foal IgG (Coming Soon)	VCF143DC	Device	5 Tests/Kit



For More Information on the
Vcheck V2400
analyzer visit:
bionote.com
customerservice@bionote.com
800-727-5169

LEARN MORE
about the V200 analyzer

